

# anti- ARMC8 antibody

#### **Product Information**

Catalog No.: FNab00587

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

None

## **Immunogen information**

Immunogen: armadillo repeat containing 8

Synonyms: None
Observed MW: 75.5kd
UniprotID: Q8IUR7

## **Application**

Reactivity: Human, Mouse, Rat Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:2000-1:10000; IHC: 1:50-1:500

Image:


#### 1

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Hela cells were subjected to SDS PAGE followed by western blot with FNab00587(ARMC8 Antibody) at dilution of 1:6000