

anti- ARP3 antibody

Product Information

Catalog No.:	FNab00598
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Background	

Background

None

Immunogen information

Immunogen:	ARP3 actin-related protein 3 homolog
Synonyms:	Actin like protein 3, Actin related protein 3, ACTR3, ACTR3B, Arp11, ARP3, Arp4
Observed MW:	47 kDa
Uniprot ID :	P61158

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:3000; IP: 1:500-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com





IP Result of anti-ARP3 (IP:FNab00598, 4ug; Detection:FNab00598 1:1000) with mouse spleen tissue lysate 4000ug.

Hela cells were subjected to SDS PAGE followed by western blot with FNab00598(ARP3 antibody) at dilution of 1:500

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889

40kd

-35kd

www.fn-test.com