

anti- ASB9 antibody

Product Information

Catalog No.: FNab00625

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Substrate-recognition component of a SCF-like ECS(Elongin-Cullin-SOCS-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins. Recognizes at least two forms of creatine kinase, CKB and CKMT1A.

Immunogen information

Immunogen: ankyrin repeat and SOCS box-containing 9

Synonyms: ASB 9, ASB9

Observed MW: 28k-32d UniprotID: Q96DX5

Application

Reactivity: Human, Mouse Tested Application: ELISA, WB

Recommended dilution: WB: 1:200-1:1000

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com

1



Image:



human testis tissue were subjected to SDS PAGE followed by western blot with FNab00625(ASB9 antibody) at dilution of 1:300