

anti- ATP12A antibody

Product Information

Catalog No.: FNab00689

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Catalyzes the hydrolysis of ATP coupled with the exchange of H(+) and K(+) ions across the plasma membrane. Responsible for potassium absorption in various tissues.

Immunogen information

Immunogen: ATPase, H+/K+ transporting, nongastric, alpha polypeptide

Synonyms: ATP12A, ATP1AL1, Proton pump

Observed MW: 100-110 kDa

Uniprot ID: P54707

Application

Reactivity: Human

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:1000

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab00689(ATP12A Antibody) at dilution of 1:500