

# anti- ATP6V0D1 antibody

#### **Product Information**

Catalog No.: FNab00714

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

Subunit of the integral membrane V0 complex of vacuolar ATPase. Vacuolar ATPase is responsible for acidifying a variety of intracellular compartments in eukaryotic cells, thus providing most of the energy required for transport processes in the vacuolar system. May play a role in coupling of proton transport and ATP hydrolysis(By similarity). May play a role in cilium biogenesis through regulation of the transport and the localization of proteins to the cilium(By similarity).

#### **Immunogen information**

Immunogen: ATPase, H+ transporting, lysosomal 38kDa, V0 subunit d1

Synonyms: ATP6D, VPATPD

Observed MW: 41 kDa UniprotID: P61421

## **Application**

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC, FC, IF, IP

## Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com


1



Recommended dilution: WB: 1:500-1:2000; IP: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100 Image:




Immunohistochemistry of paraffin-embedded human kidney using FNab00714(ATP6V0D1 antibody) at dilution of 1:100



IP Result of anti-ATP6V0D1 (IP:FNab00714, 4ug; Detection:FNab00714 1:500) with mouse testis tissue lysate 4800ug.





mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab00714(ATP6V0D1 antibody) at dilution of 1:800

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com