

anti- AUTS2 antibody

Product Information

Catalog No.:	FNab00750
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Background	

Background

None

Immunogen information

Immunogen:	autism susceptibility candidate 2
Synonyms:	KIAA0442
Observed MW:	70 kDa
UniprotID :	Q8WXX7

Application


Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB, IHC		
Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200			
Image:			

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





-40kd

-35kd

-25kd

Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab00750(AUTS2 Antibody) at dilution of 1:50

C6 cells were subjected to SDS PAGE followed by western blot with FNab00750(AUTS2 Antibody) at dilution of 1:600

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>