

anti- AXIN2 antibody

Product Information

Catalog No.: FNab00753

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Inhibitor of the Wnt signaling pathway. Down-regulates beta-catenin. Probably facilitate the phosphorylation of beta-catenin and APC by GSK3B(By similarity).

Immunogen information

Immunogen: axin 2

Synonyms: AXIL, axin 2, Axin like protein, AXIN2, Axis inhibition protein 2,

Conductin

Observed MW: 93-100 kDa Uniprot ID: Q9Y2T1

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:2000

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Rat brain were subjected to SDS PAGE followed by western blot with FNab00753(AXIN2 Antibody) at dilution of 1:1000