

anti- BAI2-Specific antibody

Product Information

Catalog No.: FNab00791

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

BAI2 belongs to the G-protein coupled receptor 2 family and LN-TM7 subfamily. BAI2 might be involved in angiogenesis inhibition. This antibody is specific to BAI2.

Immunogen information

Immunogen: brain-specific angiogenesis inhibitor 2

Synonyms: None
Observed MW: None
UniprotID: O60241

Application

Reactivity: Human

Tested Application: ELISA, IHC

Recommended dilution: IHC: 1:20-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab00791(BAI2-Specific Antibody) at dilution of 1:50