

anti-BAIAP3 antibody

Product Information

Catalog No.: FNab00795

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen: BAI1-associated protein 3

Synonyms: BAI1 associated protein 3, BAIAP3, BAP3, KIAA0734

Observed MW: 130 kDa UniprotID: 094812

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IF

Recommended dilution: WB: 1:200-1:1000; IF: 1:10-1:100

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunofluorescent analysis of SH-SY5Y cells using FNab00795(BAIAP3 antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated Goat Anti-Rabbit IgG(H+L)



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00795(BAIAP3 Antibody) at dilution of 1:300