

anti- BANP antibody

Product Information

Catalog No.: FNab00799

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Controls V(D)J recombination during T-cell development by repressing T-cell receptor(TCR) beta enhancer function. Binds to scaffold/matrix attachment region beta(S/MARbeta), an ATC-rich DNA sequence located upstream of the TCR beta enhancer. Represses cyclin D1 transcription by recruiting HDAC1 to its promoter, thereby diminishing H3K9ac, H3S10ph and H4K8ac levels. Promotes TP53 'Ser-15' phosphorylation and nuclear accumulation, which causes cell cycle arrest(By similarity).

Immunogen information

Immunogen: BTG3 associated nuclear protein

Synonyms: BANP, BEND1, Protein BANP, SMAR1, SMARBP1

Observed MW: 56kd

Uniprot ID: Q8N9N5

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, IHC, WB

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)


Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1




Recommended dilution: WB: 1:500-1:1000; IHC: 1:50-1:500

Image:



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using FNab00799(BANP antibody) at dilution of 1:200



human kidney tissue were subjected to SDS PAGE followed by western blot with FNab00799(BANP antibody) at dilution of 1:500