

anti- BBS8 antibody

Product Information

Catalog No.: FNab00824

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen: tetratricopeptide repeat domain 8

Synonyms: None
Observed MW: 50 kDa
UniprotID: Q8TAM2

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IP

Recommended dilution: WB: 1:200-1:1000

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





IP Result of anti-BBS8 (IP:FNab00824, 3ug; Detection:FNab00824 1:400) with HeLa cells lysate 2000ug.



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab00824(BBS8 antibody) at dilution of 1:300