

# anti- BCAS2 antibody

## **Product Information**

Catalog No.:	FNab00828
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### **Immunogen information**

Immunogen:	breast carcinoma amplified sequence 2
Synonyms:	DAM1;SPF27;Snt309
Observed MW:	25 kDa
Uniprot ID :	O75934

## Application


Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB, IHC		
Recommended dilution: WB: 1:500 - 1:2000; IHC: 1:50 - 1:200			
Image:			

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

1





Immunohistochemistry of paraffin-embedded human esophagus using FNab00828( BCAS2 Antibody) at dilution of 1:100 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab00828(BCAS2 antibody) at dilution of 1:1000

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)