

anti- BCL2L10 antibody

Product Information

Catalog No.: FNab00844

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Promotes cell survival. Suppresses apoptosis induced by BAX but not BAK.

Immunogen information

Immunogen: BCL2-like 10(apoptosis facilitator)

Synonyms: BCLB
Observed MW: 22 kDa
UniprotID: Q9HD36

Application

Reactivity: Human

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human lung using FNab00844(BCL2L10 antibody) at dilution of 1:50



A2780 cells were subjected to SDS PAGE followed by western blot with FNab00844(BCL2L10 antibody) at dilution of 1:200