

anti- BCL3 antibody

Product Information

Catalog No.:	FNab00848
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Contributes to the regulation of transcriptional activation of NF-kappa-B target genes. In the cytoplasm, inhibits the nuclear translocation of the NF-kappa-B p50 subunit. In the nucleus, acts as transcriptional activator that promotes transcription of NF-kappa-B target genes. Contributes to the regulation of cell proliferation(By similarity).

Immunogen information

Immunogen:	B-cell CLL/lymphoma 3
Synonyms:	BCL4, D19S37
Observed MW:	48-60 kDa
Uniprot ID :	P20749

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution	:WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com



Image:



Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using FNab00848(BCL3 Antibody) at dilution of 1:50

HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab00848(BCL3 Antibody) at dilution of 1:600

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>