

anti- Bestrophin-3 antibody

Product Information

Catalog No.:	FNab00867
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	bestrophin 3
Synonyms:	None
Observed MW:	76-78 kDa, 50-55 kDa
Uniprot ID :	Q8N1M1

Application


Reactivity:	Human, Mouse			
Tested Application:	ELISA, WB, IHC			
Recommended dilution: WB: 1:500-1:1000; IHC: 1:50-1:500				
Image:				

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com





Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using FNab00867(Bestrophin-3 antibody at dilution of 1:200



mouse colon tissue were subjected to SDS PAGE followed by western blot with FNab00867(Bestrophin-3 antibody) at dilution of 1:800

2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)