

# anti- Beta sarcoglycan antibody

#### **Product Information**

Catalog No.: FNab00878

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

### **Immunogen information**

Immunogen: sarcoglycan, beta(43kDa dystrophin-associated glycoprotein)

Synonyms: 43DAG, A3b, Beta sarcoglycan, Beta SG, LGMD2E, SGC, SGCB

Observed MW: 48kd

Uniprot ID: Q16585

## **Application**

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:1000


Image:

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





MCF7 cells were subjected to SDS PAGE followed by western blot with FNab00878(Beta sarcoglycan antibody) at dilution of 1:600

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>