

anti- BEX2 antibody

Product Information

Catalog No.:	FNab00884
Size:	100µ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^\circ\!\!C$ for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Regulator of mitochondrial apoptosis and G1 cell cycle in breast cancer. Protects the breast cancer cells against mitochondrial apoptosis and this effect is mediated through the modulation of BCL2 protein family, which involves the positive regulation of anti-apoptotic member BCL2 and the negative regulation of pro-apoptotic members BAD, BAK1 and PUMA. Required for the normal cell cycle progression during G1 in breast cancer cells through the regulation of CCND1 and CDKN1A. Regulates the level of PP2A regulatory subunit B and PP2A phosphatase activity. The homology between human BEX1 and BEX2 is 94%. This antibody detects both BEX1 and BEX2.

Immunogen information

Immunogen:	brain expressed X-linked 2
Synonyms:	BEX1, BEX1/2, BEX2, brain expressed X linked 2, hBex2, Protein BEX2
Calculated MW:	10-15 kDa
Uniprot ID :	Q9BXY8
Application	
Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com

1



Recommended dilution: WB: 1:500-1:1000; IHC: 1:20-1:200 Image:



Immunohistochemistry of paraffin-embedded human breast cancer using FNab00884(BEX2 antibody) at dilution of 1:100

	—25kd
-	——15kd
	——10kd

human brain tissue were subjected to SDS PAGE followed by western blot with FNab00884(BEX2 antibody) at dilution of 1:500

2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)