

# anti- BLVRB antibody

### **Product Information**

Catalog No.:	FNab00911
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

The final step in heme metabolism in mammals is catalyzed by the cytosolic biliverdin reductase enzymes A and B (EC 1.3.1.24).

## Immunogen information

Immunogen:	biliverdin reductase B (flavin reductase (NADPH))
Synonyms:	FLR
Observed MW:	22 kDa
Uniprot ID :	P30043

### Application

Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB, IHC		
Recommended dilution: WB: 1:500 - 1:2000; IHC: 1:50 - 1:100			
Image:			


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889


www.fn-test.com





Immunohistochemistry of paraffin-embedded human stomach using FNab00911( BLVRB Antibody) at dilution of 1:50

HepG2 cells were subjected to SDS PAGE followed by western blot with FNab00911(BLVRB antibody) at dilution of 1:1000



2

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com