

# anti- BMP4 antibody

#### **Product Information**

Catalog No.: FNab00919

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

# **Background**

Induces cartilage and bone formation. Also act in mesoderm induction, tooth development, limb formation and fracture repair. Acts in concert with PTHLH/PTHRP to stimulate ductal outgrowth during embryonic mammary development and to inhibit hair follicle induction(By similarity).

### Immunogen information

Immunogen: bone morphogenetic protein 4

Synonyms: BMP2B, DVR4

Observed MW: 50 kDa UniprotID: P12644

# **Application**

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, IHC, WB

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200

Image:


#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com

1





Immunohistochemistry of paraffin-embedded human lung using FNab00919(BMP4 antibody) at dilution of 1:50



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab00919(BMP4 Antibody) at dilution of 1:300