

anti- BRD9 antibody

Product Information

Catalog No.:	FNab00951
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Immunogen information

Immunogen:	bromodomain containing 9
Synonyms:	None
Observed MW:	67 kDa
Uniprot ID :	Q9H8M2

Application

Reactivity:	Human, Mouse			
Tested Application:	ELISA, WB, IHC			
Recommended dilution: WB: 1:500 - 1:2000; IHC: 1:50 - 1:200				
Image:				


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel :(0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com





Immunohistochemistry of paraffin-embedded human liver using FNab00951(BRD9 antibody) at dilution of 1:200



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00951(BRD9 Antibody) at dilution of 1:1000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889