

# anti- BUD13 antibody

### **Product Information**

Catalog No.:	FNab00990
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Rackground	

#### Background

None

## **Immunogen information**

Immunogen:	BUD13 homolog(S. cerevisiae)
Synonyms:	None
Observed MW:	70-75 kDa
UniprotID :	Q9BRD0

## Application


Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution	WB: 1:500-1:2000
Image:	

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





HepG2 cells were subjected to SDS PAGE followed by western blot with FNab00990( BUD13 Antibody) at dilution of 1:600

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>