

anti- C11orf57 antibody

Product Information

Catalog No.:	FNab01006
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Rackground	

Background

None

Immunogen information

Immunogen:	chromosome 11 open reading frame 57
Synonyms:	None
Observed MW:	30-34 kDa
UniprotID :	Q6ZUT1

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution	WB: 1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab01006(C11orf57 Antibody) at dilution of 1:300

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com