

anti- C15orf58 antibody

Product Information

Catalog No.:	FNab01031
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	chromosome 15 open reading frame 58
Synonyms:	None
Observed MW:	45kd
UniprotID :	Q6ZNW5

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:WB:	1:50-1:500; IF: 1:50-1:500
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab01031(C15orf58 Antibody) at dilution of 1:100