

## anti- C1orf103 antibody

### Product Information

Catalog No.:	FNab01057
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Represses the ligand-induced transcriptional activity of retinoic acid receptor alpha(RARA). This repression may occur through direct recruitment of histone deacetylases.

### Immunogen information

Immunogen:	chromosome 1 open reading frame 103
Synonyms:	FLJ11269, LRIF1, Receptor interacting factor 1, RIF1
Observed MW:	85 kDa
Uniprot ID :	Q5T3J3

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



L02 cells were subjected to SDS PAGE followed by western blot with FNaB01057(C1orf103 Antibody) at dilution of 1:1000