

## anti- C1orf26 antibody

### Product Information

Catalog No.:	FNab01063
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	chromosome 1 open reading frame 26
Synonyms:	None
Observed MW:	100-105 kDa
UniprotID :	Q5T5J6

### Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, IP
Recommended dilution:IP:	1:500-1:1000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



IP Result of anti-C1orf26 (IP:FNab01063, 5ug;  
Detection:FNab01063 1:500) with mouse testis  
tissue lysate 8000ug.