
anti- C3 / C3b / C3c antibody

Product Information

Catalog No.:	FNab01103
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	8G4
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	complement component 3
Synonyms:	ARMD9, ASP, C3, C3b, C3c, Complement C3, complement component 3, CPAMD1
Observed MW:	115 kDa
Uniprot ID :	P01024

Application


Reactivity:	Human
Tested Application:	ELISA, IHC, WB, IF
Recommended dilution:	WB: 1:1000-1:4000; IHC: 1:20-1:200; IF: 1:50-1:500
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human liver tissue slide using FNab01103(C3 Antibody) at dilution of 1:50



human plasma tissue were subjected to SDS PAGE followed by western blot with FNab01103(C3 Antibody) at dilution of 1:2000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com