

anti- C4(alpha chain) antibody

Product Information

Catalog No.:	FNab01111
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	complement component 4A(Rodgers blood group)
Synonyms:	None
Observed MW:	95 kDa
Uniprot ID :	POC0L4

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:6000; IHC: 1:50-1:500
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human liver tissue slide using FNab01111(C4 (alpha chain) antibody at dilution of 1:200



human blood tissue were subjected to SDS PAGE followed by western blot with FNab01111(C4 (alpha chain) antibody at dilution of 1:3000