

## anti- CAB39L antibody

### Product Information

Catalog No.:	FNab01168
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Component of a complex that binds and activates STK11/LKB1. In the complex, required to stabilize the interaction between CAB39/MO25(CAB39/MO25alpha or CAB39L/MO25beta) and STK11/LKB1(By similarity).

### Immunogen information

Immunogen:	calcium binding protein 39-like
Synonyms:	None
Observed MW:	39 kDa
UniprotID :	Q9H9S4

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



K-562 cells were subjected to SDS PAGE followed by western blot with FNab01168(CAB39L antibody) at dilution of 1:400