

## anti- CARNS1 antibody

### Product Information

Catalog No.:	FNab01276
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Catalyzes the synthesis of carnosine and homocarnosine. Carnosine is synthesized more efficiently than homocarnosine.

### Immunogen information

Immunogen:	ATP-grasp domain containing 1
Synonyms:	ATPGD1, KIAA1394
Observed MW:	88 kDa
UniprotID :	A5YM72

### Application


Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



HeLa cells were subjected to SDS PAGE followed by western blot with FNab01276( CARNS1 antibody) at dilution of 1:300