

anti- CCDC28B antibody

Product Information

Catalog No.:	FNab01354
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Involved in ciliogenesis. Regulates cilia length through its interaction with MAPKAP1/SIN1 but independently of mTORC2 complex. Modulates mTORC2 complex assembly and function, possibly enhances AKT1 phosphorylation. Does not seem to modulate assembly and function of mTORC1 complex.

Immunogen information

Immunogen:	coiled-coil domain containing 28B
Synonyms:	None
Observed MW:	30-35 kDa
UniprotID :	Q9BUN5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

Image:



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab01354(CCDC28B antibody) at dilution of 1:300

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com