
anti- CD100 antibody

Product Information

Catalog No.:	FNab01407
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	sema domain, immunoglobulin domain(Ig), transmembrane domain(TM) and short cytoplasmic domain,(semaphorin) 4D
Synonyms:	A8, BB18, C9orf164, CD100, coll 4, GR3, M sema G, SEMA4D, SEMAJ, Semaphorin 4D
Observed MW:	150 kDa
Uniprot ID :	Q92854

Application

Reactivity:	Human
Tested Application:	ELISA

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com