

anti- CD20 antibody

Product Information

Catalog No.:	FNab01441
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3C0
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	membrane-spanning 4-domains, subfamily A, member 1
Synonyms:	B lymphocyte antigen CD20, B1, Bp35, CD20, LEU 16, MS4A1, MS4A1,CD20, MS4A2, S7
Observed MW:	33 kDa
Uniprot ID :	P11836

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, FC
Recommended dilution:WB:	1:2000-1:15000; IHC: 1:1000-1:15000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab01441(CD20 Antibody) at dilution of 1:10000



Raji cells were subjected to SDS PAGE followed by western blot with FNab01441(CD20 Antibody) at dilution of 1:8000