

anti- CD31 antibody

Product Information

Catalog No.:	FNab01464
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2D10
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	platelet/endothelial cell adhesion molecule
Synonyms:	CD31, EndoCAM, GPIIA, PECA1, PECAM 1, PECAM1, VEC marker
Observed MW:	130 kDa
Uniprot ID :	P16284

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:2000-1:6000; IHC: 1:500-1:3000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab01464(CD31 Antibody) at dilution of 1:500



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab01464(CD31 antibody) at dilution of 1:1000