

anti-Coilin antibody

Product Information

Catalog No.:	FNab01833
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	coilin
Synonyms:	None
Observed MW:	80 kDa
UniprotID :	P38432

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, FC, IP, IF
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275


Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab01833(COIL antibody) at dilution of 1:100



IP Result of anti-Coilin (IP:FNab01833, 3ug; Detection:FNab01833 1:500) with HEK-293 cells lysate 1000ug.



human brain tissue were subjected to SDS PAGE followed by western blot with FNab01833(COIL antibody) at dilution of 1:800

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com