

anti- CPI17 antibody

Product Information

Catalog No.:	FNab01919
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	protein phosphatase 1, regulatory(inhibitor) subunit 14A
Synonyms:	None
Observed MW:	17 kDa
UniprotID :	Q96A00

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:1000-1:4000; IF: 1:10-1:100
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

HepG2 cells were subjected to SDS PAGE followed by western blot with FNab01919(CPI17 antibody) at dilution of 1:1000

