

## anti- CPLX3 antibody

### Product Information

Catalog No.:	FNab01923
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Positively regulates a late step in synaptic vesicle exocytosis.

### Immunogen information

Immunogen:	complexin 3
Synonyms:	None
Observed MW:	20-23 kDa
UniprotID :	Q8WVH0

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	



### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)

IP Result of anti-CPLX3 (IP:FNab01923, 3ug;  
Detection:FNab01923 1:300) with mouse brain  
tissue lysate 6000ug.



mouse brain tissue were subjected to SDS PAGE  
followed by western blot with FNab01923(CPLX3  
antibody) at dilution of 1:300