

anti- CSGALNACT2 antibody

Product Information

Catalog No.:	FNab02014
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Transfers 1,4-N-acetylgalactosamine(GalNAc) from UDP-GalNAc to the non-reducing end of glucuronic acid(GlcUA). Required for addition of the first GalNAc to the core tetrasaccharide linker and for elongation of chondroitin chains.

Immunogen information

Immunogen:	chondroitin sulfate N-acetylgalactosaminyltransferase 2
Synonyms:	CHGN2, GALNACT2
Observed MW:	63-70 kDa
UniprotID :	Q8N6G5

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



K-562 cells were subjected to SDS PAGE followed by western blot with FNab02014(CSGALNACT2 antibody) at dilution of 1:800