

anti- CXCL17 antibody

Product Information

Catalog No.:	FNab02097
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Plays a role in angiogenesis and possibly in the development of tumors. May be a housekeeping chemokine regulating recruitment of nonactivated blood monocytes and immature dendritic cells into tissues. May play a role in the innate defense against infections.

Immunogen information

Immunogen:	chemokine(C-X-C motif) ligand 17
Synonyms:	C X C motif chemokine 17, CXCL17, Dcip1, DMC, UNQ473, VCC 1, VCC1, VEGF co regulated chemokine 1
Observed MW:	None
Uniprot ID :	Q6UXB2

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:20-1:200

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com

Image:



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab02097(CXCL17 Antibody) at dilution of 1:50