

anti- Cyclin A1 antibody

Product Information

Catalog No.:	FNab02120
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May be involved in the control of the cell cycle at the G1/S(start) and G2/M(mitosis) transitions. May primarily function in the control of the germline meiotic cell cycle and additionally in the control of mitotic cell cycle in some somatic cells.

Immunogen information

Immunogen:	cyclin A1
Synonyms:	CCNA1, cyclin A1
Observed MW:	65 kDa
Uniprot ID :	P78396

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab02120(CCNA1 Antibody) at dilution of 1:50



SKOV-3 cells were subjected to SDS PAGE followed by western blot with FNab02120(CCNA1 antibody) at dilution of 1:500