

---

## Cyclin D3 antibody

### Product Information


Catalog No.:	FNab02126
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	cyclin D3
Synonyms:	CCND3, cyclin D3, G1/S-specific cyclin-D3
Observed MW:	33 kDa
Uniprot ID :	P30281

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using FNab02126(Cyclin D3 Antibody) at dilution of 1:200



C6 cells were subjected to SDS PAGE followed by western blot with FNab02126(Cyclin D3 Antibody) at dilution of 1:1000