
anti- CYP2C8 antibody

Product Information

Catalog No.:	FNab02156
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	cytochrome P450, family 2, subfamily C, polypeptide 8
Synonyms:	CPC8, CYP2C8, CYP11C8, Cytochrome P450 2C8
Uniprot ID :	P10632

Application


Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:500-1:1000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human liver cancer using FNab02156(CYP2C8 antibody) at dilution of 1:500