

## anti- DAB2 antibody

### Product Information

Catalog No.:	FNab02228
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	disabled homolog 2, mitogen-responsive phosphoprotein(Drosophila)
Synonyms:	DAB2, Disabled homolog 2, DOC 2, DOC2, FLJ26626
Observed MW:	96 kDa
Uniprot ID :	P98082

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF, FC, IP
Recommended dilution:	WB: 1:500-1:5000; IP: 1:500-1:2000; IF: 1:20-1:200; IHC: 1:50-1:200
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889


[www.fn-test.com](http://www.fn-test.com)




Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab02228(DAB2 Antibody) at dilution of 1:100



Immunofluorescent analysis of HepG2 cells, using DAB2 antibody FNab02228 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



IP Result of anti-DAB2 (IP:FNab02228, 4ug;  
Detection:FNab02228 1:2000) with HeLa cells  
lysate 4000ug.



HeLa cells were subjected to SDS PAGE followed  
by western blot with FNab02228(DAB2 antibody)  
at dilution of 1:1000