

anti- DAZL antibody

Product Information

Catalog No.:	FNab02252
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

RNA-binding protein, which is essential for gametogenesis in both males and females. Plays a central role during spermatogenesis. Acts by binding to the 3'-UTR of mRNA, specifically recognizing GUU triplets, and thereby regulating the translation of key transcripts (By similarity).

Immunogen information

Immunogen:	deleted in azoospermia-like
Synonyms:	DAZH, DAZL1, DAZLA, SPGYLA
Observed MW:	35kd
UniprotID :	Q92904

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP Result of anti-DAZL (IP:FNab02252, 4ug;
Detection:FNab02252 1:500) with mouse testis
tissue lysate 4400ug.



mouse testis tissue were subjected to SDS PAGE
followed by western blot with FNab02252(DAZL
antibody) at dilution of 1:400