

anti-Desmin antibody

Product Information

Catalog No.:	FNab02342
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	6H3
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	desmin
Synonyms:	CMD1I, CSM1, CSM2, DES, desmin, FLJ12025, FLJ39719, FLJ41013, FLJ41793
Observed MW:	50-55 kDa
UniprotID :	P17661

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human heart using FNab02342(DES antibody) at dilution of 1:50



human heart tissue were subjected to SDS PAGE followed by western blot with FNab02342(DES antibody) at dilution of 1:1000