

anti- DGCR2 antibody

Product Information

Catalog No.:	FNab02350
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Putative adhesion receptor, that could be involved in cell-cell or cell-matrix interactions required for normal cell differentiation and migration.

Immunogen information

Immunogen:	DiGeorge syndrome critical region gene 2
Synonyms:	IDD, KIAA0163
Observed MW:	61-65 kDa
Uniprot ID :	P98153

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:50-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com