

anti- DNAJB2 antibody

Product Information

Catalog No.:	FNab02448
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2H8
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	DnaJ(Hsp40) homolog, subfamily B, member 2
Synonyms:	HSJ1, HSPF3
Observed MW:	38 kDa
Uniprot ID :	P25686

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:1000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



K-562 cells were subjected to SDS PAGE followed by western blot with FNab02448(DNAJB2 antibody) at dilution of 1:300