

## anti- DRG1 antibody

### Product Information

|               |                                                                                                                 |
|---------------|-----------------------------------------------------------------------------------------------------------------|
| Catalog No.:  | FNab02537                                                                                                       |
| Size:         | 100µg                                                                                                           |
| Form:         | liquid                                                                                                          |
| Purification: | Immunogen affinity purified                                                                                     |
| Purity:       | ≥95% as determined by SDS-PAGE                                                                                  |
| Host:         | Rabbit                                                                                                          |
| Clonality:    | polyclonal                                                                                                      |
| Clone ID:     | None                                                                                                            |
| IsoType:      | IgG                                                                                                             |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Background

Critical regulator of cell growth under specific conditions. Implicated in differentiation and cell cycle arrest.

### Immunogen information

|              |                                                 |
|--------------|-------------------------------------------------|
| Immunogen:   | developmentally regulated GTP binding protein 1 |
| Synonyms:    | NEDD3                                           |
| Observed MW: | 41 kDa                                          |
| UniprotID :  | Q9Y295                                          |

### Application


|                       |                                  |
|-----------------------|----------------------------------|
| Reactivity:           | Human, Mouse, Rat                |
| Tested Application:   | ELISA, WB, IF                    |
| Recommended dilution: | WB: 1:500-1:2000; IF: 1:20-1:200 |
| Image:                |                                  |

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



K-562 cells were subjected to SDS PAGE followed by western blot with FNab02537(DRG1 antibody) at dilution of 1:4000