

## anti- DSG4 antibody

### Product Information

Catalog No.:	FNab02545
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Component of intercellular desmosome junctions. Involved in the interaction of plaque proteins and intermediate filaments mediating cell-cell adhesion. Coordinates the transition from proliferation to differentiation in hair follicle keratinocytes (By similarity).

### Immunogen information

Immunogen:	desmoglein 4
Synonyms:	CDHF13
Observed MW:	160 kDa
UniprotID :	Q86SJ6

### Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



mouse skin tissue were subjected to SDS PAGE followed by western blot with FNab02545( DSG4 Antibody) at dilution of 1:600