

## anti- DTD1 antibody

### Product Information

Catalog No.:	FNab02548
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

ATPase involved in DNA replication, may facilitate loading of CDC45 onto pre-replication complexes. May hydrolyze D-tyrosyl-tRNA(Tyr) into D-tyrosine and free tRNA(Tyr), a possible defense mechanism against a harmful effect of D-tyrosine.

### Immunogen information

Immunogen:	D-tyrosyl-tRNAdeacylase 1 homolog(S. cerevisiae)
Synonyms:	C20orf88, DUEB, HARS2
Observed MW:	23-28kd
UniprotID :	Q8TEA8

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



human brain tissue were subjected to SDS PAGE followed by western blot with FNab02548(DTD1 antibody) at dilution of 1:500