

## anti- DTWD2 antibody

### Product Information

|               |   |
|---------------|---|
| Catalog No.:  | FNab02552   |
| Size:         | 100µg   |
| Form:         | liquid  |
| Purification: | Immunogen affinity purified   |
| Purity:       | ≥95% as determined by SDS-PAGE  |
| Host:         | Rabbit  |
| Clonality:    | polyclonal  |
| Clone ID:     | None  |
| IsoType:      | IgG   |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Background

None

### Immunogen information

|              |                         |
|--------------|-------------------------|
| Immunogen:   | DTW domain containing 2 |
| Synonyms:    | None                    |
| Observed MW: | 33 kDa                  |
| UniprotID :  | Q8NBA8                  |

### Application


|                          |              |
|--------------------------|--------------|
| Reactivity:              | Human        |
| Tested Application:      | ELISA, IP    |
| Recommended dilution:IP: | 1:200-1:2000 |
| Image:                   |              |

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



IP Result of anti-DTWD2 (IP:FNab02552, 3ug;  
Detection:FNab02552 1:300) with COLO 320 cells  
lysate 2000ug.